

The Synd-EZ kit provides a solution for damage to the syndesmosis joint. This simple yet effective device is designed to stabilize the ankle and can be used with or without a fracture plate.

Versatile implant design


- Ideal for use with or without a fracture plate
- Washer included for use without a fracture plate*
- Top hat is low profile and self centering on fracture plates

Adjustable knotless construct

- Simplifies surgical procedure
- Offers confidence of flipped medial button from white/black toggle suture

Counter traction suture & tension handle

- Integrated for predictably smooth tension of construct




* Washer is only available in the Synd-EZ Titanium Kit 11223.

Synd-EZ Kit *Syndesmosis Repair*

SYND-EZ KIT COMPONENTS

1. Top hat
2. White/blue adjustment suture
3. Solid blue counter-traction suture
4. Button
5. Solid white passing suture
6. White/black toggle suture
7. 3.5mm cannulated drill bit
8. 2.4mm passing pin
9. 1.5mm K-wire
10. Black tension handle
11. Washer*



SYND-EZ ORDERING INFORMATION

Synd-EZ			
Part #	Width of Button (mm)	Length of Button (mm)	Description
11223	3.2	12	Synd-EZ Titanium Kit (Ti)
11224	3.2	12	Synd-EZ Stainless Steel Kit (SS)

Synd-EZ Instrumentation (for use with a fracture plate)			
Part #	Description	Sterile/Non-sterile	Single-use/Reusable
11272	Synd-EZ Drill Guide	Non-sterile	Reusable

* Washer is only available in the Synd-EZ Titanium Kit 11223.

Anika Therapeutics, Inc.

Parcus Medical, LLC
 6423 Parkland Dr, Sarasota, FL 34243
 1-877-746-2972 | ParcusCS@anika.com
www.anika.com | Anika. Restore Active Living.™ | Stay Active®

This document and information is intended for markets where regulatory approval has been granted.
 Anika, Parcus, Synd-EZ, Stay Active, and Restore Active Living are trademarks and/or registered trademarks of Anika Therapeutics, Inc. and its affiliates in certain jurisdictions.
 System designed and manufactured in the USA | Printed in the USA
 ©2023 Anika Therapeutics, Inc. All rights reserved.
 QD 8042 | AML-900-227 REV 01