

Parcus Braid® is available in multiple configurations: suture, suture tape, and infinity loop. It is recommended for use in soft tissue approximation and allograft tissues in sports medicine procedures.

Features & Benefits

Made of Ultra High Molecular Weight Polyethylene (UHMWPE)


- High strength
- Non-absorbable

Variety of color combinations

- Suture management in complicated repairs

Various needle options

- Ideal for open surgical procedures



Parcus Braid Suture					
Part #	Size	Description	Total Length (in)	Working/Loop Length (in)	QTY/ Box
20046	#1	Parcus Braid suture (wht/blu)	36	-	12
20050	#2	Parcus Braid suture (wht/blu)	36	-	12
20049	#2	Two Parcus Braid sutures (wht/blu & wht/blk)	36	-	12
20049S	#2	Two Parcus Braid sutures (wht/blu & wht/blk)	36	-	1
20035S	#2	Parcus Braid suture (blu), w/needles (MO-6 & CE)	36	-	1
20055	#2	Parcus Braid suture (wht), w/needles (MO-6)	36	-	12
20055S	#2	Parcus Braid suture (wht), w/needles (MO-6)	36	-	1
20054	#2	Parcus Braid suture (wht/blk), w/needles (MO-6)	36	-	12
20054S	#2	Parcus Braid suture (wht/blk), w/needles (MO-6)	36	-	1
20059	#5	Parcus Braid suture (wht), w/needles (CCS)	36	-	12
20059S	#5	Parcus Braid suture (wht), w/needles (CCS)	36	-	1

Parcus Braid Suture Tape					
Part #	Size	Description	Total Length (in)	Working/Loop Length (in)	QTY/ Box
20032S	1.6mm	Parcus Braid suture tape (wht/blu)	39	27	1
20034S	1.6mm	Parcus Braid suture tape (wht/blk)	39	27	1
20042	2mm	Parcus Braid suture tape (wht/blu)	39	7	6
20042S	2mm	Parcus Braid suture tape (wht/blu)	39	7	1

Parcus Braid Infinity Loop					
Part #	Size	Description	Total Length (in)	Working/Loop Length (in)	QTY/ Box
20053	#2	Parcus Braid infinity loop (wht), w/needle (ST)	24	20	12
20053S	#2	Parcus Braid infinity loop (wht), w/needle (ST)	24	20	1
20056	#2	Parcus Braid infinity loop (wht/blu), w/needle (ST)	24	20	12
20056S	#2	Parcus Braid infinity loop (wht/blu), w/needle (ST)	24	20	1
20060	#5	Parcus Braid infinity loop (wht/blu), w/needle (ST)	24	20	12
20188S	#5	Parcus Braid infinity loop (wht/blu), w/needle (CCS)	24	20	1
20174S	1.6mm	Parcus Braid suture tape, infinity loop (wht/blu), w/needle (ST)	24	20	1

Key: Needles	
Type	Description
CE	36mm, curved, reverse cutting
MO-6	26mm, curved, taper point
MO-6 & CE	MO-6 needle on one side, CE needle on the other side
BB-1	17mm, curved, taper point
CCS	48mm, curved, conventional cutting
ST	60mm, straight, taper point