

# Twist & Triple Twist

## PEEK Screw-in Suture Anchors

Twist & Triple Twist PEEK Screw-In Suture Anchors are recommended for use in both large and small-joint repairs.

### Features & Benefits

#### Fully-threaded design

- Optimal cortical fixation
- Reduces the risk of "pull-back"

#### Made of PEEK

- Non-absorbable, radiolucent, MR safe



Twist PEEK Screw-In Suture Anchors					
Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10511T	4.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (blu & wht/grn)	PEEK	1
10517T	5.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (blu & wht/grn)	PEEK	1
10518T	6.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (blu & wht/grn)	PEEK	1
10521	4.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (wht & wht/blu), w/needles (MO-6)	PEEK	1
10522	5.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (wht & wht/blu), w/needles (MO-6)	PEEK	1
10523	6.5	12.6	Twist, screw-in, suture anchor, w/2, #2 sutures (wht & wht/blu), w/needles (MO-6)	PEEK	1

Triple Twist PEEK Screw-In Suture Anchors					
Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10524T	5.5	12.6	Triple Twist, screw-in, suture anchor, w/3, #2 sutures (blu, wht/grn & wht/blk)	PEEK	1
10525T	6.5	12.6	Triple Twist, screw-in, suture anchor, w/3, #2 sutures (blu, wht/grn & wht/blk)	PEEK	1

Twist PEEK Screw-In Suture Anchors w/Two Suture Tapes					
Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11003	4.5	12.6	Twist, screw-in, suture anchor, w/2 1.6mm suture tapes (wht/blu & wht/blk)	PEEK	1
11283	4.5	12.6	Twist, screw-in, suture anchor, w/2 1.6mm suture tapes (blu, wht/blu)	PEEK	1
11285	4.5	12.6	Twist, screw-in, suture anchor, w/2 1.6mm suture tapes (blk, wht/blk)	PEEK	1
11004	5.5	12.6	Twist, screw-in, suture anchor, w/2 1.6mm suture tapes (wht/blu & wht/blk)	PEEK	1
11284	5.5	12.6	Twist, screw-in, suture anchor, w/2 1.6mm suture tapes (blu, wht/blu)	PEEK	1
11286	5.5	12.6	Twist, screw-in, suture anchor, w/2 1.6mm suture tapes (blk, wht/blk)	PEEK	1
11005	6.5	12.6	Twist, screw-in, suture anchor, w/2 1.6mm suture tapes (wht/blu & wht/blk)	PEEK	1

Twist PEEK SST (Static Suture Tape) Screw-In Suture Anchors w/One Suture Tape					
Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11349	4.5	12.6	Twist SST, screw-in, suture anchor, w/1, 2mm suture tape (blu)	PEEK	1
11350	4.5	12.6	Twist SST, screw-in, suture anchor, w/1, 2mm suture tape (blk)	PEEK	1
11351	5.5	12.6	Twist SST, screw-in, suture anchor, w/1, 2mm suture tape (blu)	PEEK	1
11352	5.5	12.6	Twist SST, screw-in, suture anchor, w/1, 2mm suture tape (blk)	PEEK	1

Instrumentation for Twist PEEK Screw-In Suture Anchors					
Part #	Diameter (mm)	Joint	Description	Sterile/Non-Sterile	Single-use/Reusable
10365	3.2	Shoulder	Awl	Non-Sterile	Reusable
10338	4.1	Shoulder	Awl	Non-Sterile	Reusable
10609	4.5	Shoulder	Punch/tap	Non-Sterile	Reusable
10277	5.5	Shoulder	5.5mm or 6.5mm punch/tap	Non-Sterile	Reusable
10222	5.5	Shoulder	5.5mm or 6.5mm punch/tap, hudson	Non-Sterile	Reusable
10505U	3.2	Small Joint	Drill bit, w/positive stop	Non-Sterile	Reusable
10540U	3.2	Shoulder	Drill bit, w/positive stop	Non-Sterile	Reusable
10541U	4	Shoulder	Drill bit	Non-Sterile	Reusable
11000U	4	Small Joint	Drill bit	Non-Sterile	Reusable
11241	4	Small Joint	4.0 x 100mm drill guide, 6-point	Non-Sterile	Reusable

### Anika Therapeutics, Inc.

Parcus Medical, LLC  
6423 Parkland Dr, Sarasota, FL 34243  
1-877-746-2972 | [ParcusCS@anika.com](mailto:ParcusCS@anika.com)

[www.anika.com](http://www.anika.com) | Anika. Restore Active Living.™ | Stay Active®