

Titanium Interference Screws

Titanium Interference Screws are recommended for the attachment of soft tissue to bone. They are indicated for use in the fixation of ligaments and tendons in patients requiring ligament or tendon repair.



Features & Benefits

Made of 6Al-4V ELI titanium alloy

- Provides the strength and stiffness required for successful ligament reconstruction

Fully-threaded design

- Optimal fixation for both soft tissue and BTB grafts

Rounded edge

- Protects soft tissue grafts

| Titanium Interference Screws | | | | | |
|------------------------------|---------------|-------------|--------------------|----------|---------|
| Part # | Diameter (mm) | Length (mm) | Description | Material | QTY/Box |
| 10001 | 7 | 20 | Interference screw | Titanium | 1 |
| 10002 | 7 | 25 | Interference screw | Titanium | 1 |
| 10003 | 7 | 30 | Interference screw | Titanium | 1 |
| 10004 | 7 | 35 | Interference screw | Titanium | 1 |
| 10005 | 8 | 20 | Interference screw | Titanium | 1 |
| 10006 | 8 | 25 | Interference screw | Titanium | 1 |
| 10007 | 8 | 30 | Interference screw | Titanium | 1 |
| 10008 | 8 | 35 | Interference screw | Titanium | 1 |
| 10009 | 9 | 20 | Interference screw | Titanium | 1 |
| 10010 | 9 | 25 | Interference screw | Titanium | 1 |
| 10011 | 9 | 30 | Interference screw | Titanium | 1 |
| 10012 | 9 | 35 | Interference screw | Titanium | 1 |
| 10014 | 10 | 25 | Interference screw | Titanium | 1 |
| 10015 | 10 | 30 | Interference screw | Titanium | 1 |
| 10016 | 10 | 35 | Interference screw | Titanium | 1 |
| 10163 | 11 | 25 | Interference screw | Titanium | 1 |
| 10164 | 11 | 30 | Interference screw | Titanium | 1 |
| 10165 | 11 | 35 | Interference screw | Titanium | 1 |
| 10018 | 12 | 25 | Interference screw | Titanium | 1 |
| 10019 | 12 | 30 | Interference screw | Titanium | 1 |
| 10020 | 12 | 35 | Interference screw | Titanium | 1 |

| Instrumentation for Titanium Interference Screws | | | | |
|--|---------------|---|----------------------|----------------------|
| Part # | Diameter (mm) | Description | Sterile/ Non-sterile | Single-use/ Reusable |
| 10838 | 4.3 | Tunnel notcher | Non-sterile | Reusable |
| 10274 | 3.5 | Hex titanium interference screw driver shaft, w/1.5mm cannulation | Non-sterile | Reusable |
| 10041 | 3.5 | Hex titanium interference screw driver shaft, hudson, w/1.5mm cannulation | Non-sterile | Reusable |
| 10372 | 6.4 | Interference screw threaded dilator, w/1.5mm cannulation | Non-sterile | Reusable |
| G900037 | - | Ratcheting handle, hudson | Non-sterile | Reusable |
| 10499 | 1.5 | Pin Pack | Sterile | Single Use |