

AcTiFlip is a button fixation system recommended for proximal biceps tenodesis and distal biceps repair. The titanium AcTiFlip button is designed for unicortical and bicortical techniques. AcTiFlip is available in two configurations. The innovative “Cinch Loop” option saves time in the operating room. The “Infinity Loop” option is designed for conventional whip stitching of the biceps tendon.

Features & Benefits

AcTiFlip Inserter

- Robust and reliable insertion of the button

Retractable finger

- Ability to reload the button if needed

Cinch loop

- Eliminates the need to whip stitch the biceps tendon
- Precise control and tension on the biceps tendon




AcTiFlip Biceps Tenodesis System

Part #	Width of Button (mm)	Length of Button (mm)	Description	Length of Loop (mm)	QTY/Box
11183SB	3	10	AcTiFlip SB, cinch loop, w/1, #2 suture (wht/blu), w/ needles (MO-6)	-	1
11185SB	3	10	AcTiFlip SB, w/1 whip (wht/blu), w/suture passer	-	1

AcTiFlip Instrumentation

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
11152SB	3.5	Shoulder	AcTiFlip SB, drill bit	Sterile	Single-use
11220	3.5	Shoulder	Cannulated reamer w/1.5mm K-Wire	Sterile	Single-use



AcTiFlip Components

1	AcTiFlip button
2	Inserter shaft with post
3	Retractable finger
4	Handle with activation button
5	Reload pin
6	AcTiFlip cinch loop
7	AcTiFlip drill bit

Parcus Medical has joined **Anika**

PARCUS

Parcus Medical, LLC.

6423 Parkland Dr, Sarasota, FL

1-941-755-7965 • www.parcusmedical.com

ParcusCS@anika.com