

SLiK is a knotless screw-in anchor recommended for use in both large and small-joint repairs. Combining the SLiK anchor with Parcus Braid Suture Tape helps to create an ideal knotless double row rotator cuff repair. Available in natural PEEK-OPTIMA.

## Features & Benefits

### Activation button

- Allows for one-handed insertion which saves time in the operating room

### Large suture guide

- Ability to set tension prior to fixation
- Accommodates multiple suture/tape configurations

### Knotless

- Provides a step-saving alternative to conventional “knotted” suture anchors
- Eliminates “knot stacks” associated with soft tissue irritation



### SLiK PEEK Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10915	4.75	15	SLiK, knotless, screw-in, suture anchor, w/suture passer	PEEK	1
10918	5.5	15	SLiK, knotless, screw-in, suture anchor, w/suture passer	PEEK	1
10919	6.25	15	SLiK, knotless, screw-in, suture anchor, w/suture passer	PEEK	1

### Instrumentation for SLiK Anchors

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
11092	4.1	Shoulder	SLiK, awl	Unsterile	Reusable
11000U	4	Small Joint	Drill bit	Unsterile	Reusable
11241	4	Small Joint	4.0 x 100mm drill guide, 6-point	Unsterile	Reusable



Parcus Medical has joined **Anika**



**Parcus Medical, LLC.**

6423 Parkland Dr, Sarasota, FL

1-941-755-7965 · [www.parcusmedical.com](http://www.parcusmedical.com)

[ParcusCS@anika.com](mailto:ParcusCS@anika.com)