

GFS Ultimate is a suspensory fixation device recommended for use in the fixation of ligaments and tendons in patients requiring ligament or tendon repairs. It combines the convenience of adjustability with the strength and stiffness of conventional fixed-length loop devices.

Features & Benefits

Adjustable loop (Min:12mm to Max:60mm)

- Simplifies surgical procedure
- Eliminates the need to stock multiple sizes
- Can be both tightened and loosened

Loop made of Ultra High Molecular Weight Polyethylene (UHMWPE)

- High strength
- Non-absorbable

Button made of 6Al-4V ELI titanium alloy

- Provides the strength and stiffness required for successful ligament reconstruction

Designed to be used with various surgical techniques/instrumentation

- Ideal for both transtibial and anterior medial approaches
- Ideal for both outside-in instrumentation and reverse cutting tunnel instrumentation
- Ideal for stepped femoral tunnels or a single diameter femoral tunnel



GFS Ultimate

Part #	Width of Button (mm)	Length of Button (mm)	Description	Length of Loop	QTY/Box
10972	3.9	12	GFS, ultimate mini	Min: 12, Max: 60	1
10974	3.9	16.5	GFS, ultimate large	Min: 12, Max: 60	1

Instrumentation for GFS Ultimate

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10057	-	Knee	GFS Depth Gauge	Unsterile	Reusable
10613U	5	Knee	Headed reamer, w/2.5mm cannulation ¹	Unsterile	Reusable

Parcus Medical has joined **Anika**

PARCUS

Parcus Medical, LLC.

6423 Parkland Dr, Sarasota, FL

1-941-755-7965 • www.parcusmedical.com

ParcusCS@anika.com

©2021 Parcus Medical, LLC. | QD 8032 Rev 3 3/2021

Parcus is a registered trademark of Parcus Medical, LLC.

¹Use w/ GFS Ultimate/GFS Ultimate Mini

