

GFS Naked is a suspensory fixation device recommended for use in the fixation of ligaments and tendons in patients requiring ligament or tendon repairs. The larger button may eliminate the stepped tunnel technique required by smaller buttons. It is an excellent choice for revision surgeries.

Features & Benefits

Button made of 6Al-4V ELI titanium alloy

- Provides the strength and stiffness required for successful ligament reconstruction

GFS Mini					
Part #	Width of Button (mm)	Length of Button (mm)	Description	Length of Loop (mm)	QTY/Box
10787	1.4	5.9	GFS, naked nano	-	1
10788	2	8	GFS, naked micro	-	1
10789	3.9	12	GFS, naked mini	-	1
10790	3.9	14	GFS, naked standard	-	1
10791	4.4	16.5	GFS, naked large	-	1
10792	4.4	21	GFS, naked XL	-	1



Parcus Medical has joined **Anika**

PARCUS

Parcus Medical, LLC.

6423 Parkland Dr, Sarasota, FL

1-941-755-7965 • www.parcusmedical.com

ParcusCS@anika.com