



Titanium Interference Screws are recommended for the attachment of soft tissue to bone.

Features & Benefits

Made of 6Al-4V ELI titanium alloy

- Provides the strength and stiffness required for successful ligament reconstruction

Fully-threaded design

- Optimal fixation for both soft tissue and BTB grafts

Rounded edge

- Protects soft tissue grafts

Titanium Interference Screws					
Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10001	7	20	Interference screw	Titanium	1
10002	7	25	Interference screw	Titanium	1
10003	7	30	Interference screw	Titanium	1
10004	7	35	Interference screw	Titanium	1
10005	8	20	Interference screw	Titanium	1
10006	8	25	Interference screw	Titanium	1
10007	8	30	Interference screw	Titanium	1
10008	8	35	Interference screw	Titanium	1
10009	9	20	Interference screw	Titanium	1
10010	9	25	Interference screw	Titanium	1
10011	9	30	Interference screw	Titanium	1
10012	9	35	Interference screw	Titanium	1
10014	10	25	Interference screw	Titanium	1
10015	10	30	Interference screw	Titanium	1
10016	10	35	Interference screw	Titanium	1
10163	11	25	Interference screw	Titanium	1
10164	11	30	Interference screw	Titanium	1
10165	11	35	Interference screw	Titanium	1
10018	12	25	Interference screw	Titanium	1
10019	12	30	Interference screw	Titanium	1
10020	12	35	Interference screw	Titanium	1

Instrumentation for Titanium Interference Screws					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10838	4.3	Knee	Tunnel notcher	Unsterile	Reusable
10274	3.5	Knee	Hex titanium interference screw driver shaft, w/1.5mm cannulation	Unsterile	Reusable
10041	3.5	Knee	Hex titanium interference screw driver shaft, hudson, w/1.5mm cannulation	Unsterile	Reusable
10372	6.4	Knee	Interference screw threaded dilator, w/1.5mm cannulation	Unsterile	Reusable



Knee