

Suture, Suture Tape & Infinity Loop

Anika offers standard sutures, suture tapes, and infinity loop sutures in multiple configurations and color options. These sutures are recommended for use in soft tissue approximation and allograft tissues in sports medicine procedures.

Features & Benefits

Made with Ultra High Molecular Weight Polyethylene (UHMWPE)


- High strength
- Non-absorbable
- New AnikaBraid™ 2mm Suture Tapes are a 50/50 blend of UHMWPE and Polyester

Variety of pattern and color combinations

- Allows suture management in complicated repairs
- Easily visible in arthroscopic setting
- Suture identifiable to surgeons with color vision deficiency

Various needle options

- Ideal for open surgical procedures



Suture				
Part #	Size	Description	Total Length (in)	QTY/ Box
20046	#1	ParcusBraid® suture (wht/blu)	36	12
20050	#2	ParcusBraid suture (wht/blu)	36	12
20049	#2	Two ParcusBraid sutures (wht/blu & wht/blk)	36	12
20049S	#2	Two ParcusBraid sutures (wht/blu & wht/blk)	36	1
20035S	#2	ParcusBraid suture (blu), w/needles (MO-6 & CE)	36	1
20055	#2	ParcusBraid suture (wht), w/needles (MO-6)	36	12
20055S	#2	ParcusBraid suture (wht), w/needles (MO-6)	36	1
20054	#2	ParcusBraid suture (wht/blk), w/needles (MO-6)	36	12
20054S	#2	ParcusBraid suture (wht/blk), w/needles (MO-6)	36	1
20059	#5	ParcusBraid suture (wht), w/needles (CCS)	36	12
20059S	#5	ParcusBraid suture (wht), w/needles (CCS)	36	1

Suture Tape					
Part #	Size	Description	Total Length (in)	Working/Loop Length (in)	QTY/ Box
20032S	1.6mm	ParcusBraid suture tape (wht/blu)	39	27	1
20034S	1.6mm	ParcusBraid suture tape (wht/blk)	39	27	1
20042	2mm	ParcusBraid suture tape (wht/blu)	40	32	6
20042S	2mm	ParcusBraid suture tape (wht/blu)	40	32	1
20238	2mm	AnikaBraid suture tape (wht/blu)	34	10	12
20238S	2mm	AnikaBraid suture tape (wht/blu)	34	10	1
20239	2mm	AnikaBraid suture tape (blu)	34	10	12
20239S	2mm	AnikaBraid suture tape (blu)	34	10	1
20240	2mm	AnikaBraid suture tape (blk/blu)	34	10	12
20240S	2mm	AnikaBraid suture tape (blk/blu)	34	10	1

Infinity Loop					
Part #	Size	Description	Total Length (in)	Working/Loop Length (in)	QTY/ Box
20053	#2	ParcusBraid infinity loop (wht), w/needle (ST)	24	20	12
20053S	#2	ParcusBraid infinity loop (wht), w/needle (ST)	24	20	1
20056	#2	ParcusBraid infinity loop (wht/blu), w/needle (ST)	24	20	12
20056S	#2	ParcusBraid infinity loop (wht/blu), w/needle (ST)	24	20	1
20060	#5	ParcusBraid infinity loop (wht/blu), w/needle (ST)	24	20	12
20188S	#5	ParcusBraid infinity loop (wht/blu), w/needle (CCS)	24	20	1
20174S	1.6mm	ParcusBraid suture tape, infinity loop (wht/blu), w/needle (ST)	24	20	1

Anika Therapeutics, Inc.

Parcus Medical, LLC

6423 Parkland Dr, Sarasota, FL 34243

1-888-721-1600 | customerservice@anika.com

www.anika.com | Anika. Restore Active Living.® | Stay Active®

This document and information is intended for markets where regulatory approval has been granted.

Anika, Parcus, AnikaBraid, ParcusBraid, Stay Active, and Restore Active Living are trademarks and/or registered trademarks of Anika Therapeutics, Inc. and its affiliates in certain jurisdictions.

System designed and manufactured in the USA | Printed in the USA

©2024 Anika Therapeutics, Inc. All rights reserved.

AML-900-364 REV 02

Key: Needles

Type	Description
CE	36mm, curved, reverse cutting
MO-6	26mm, curved, taper point
MO-6 & CE	MO-6 needle on one side, CE needle on the other side
BB-1	17mm, curved, taper point
CCS	48mm, curved, conventional cutting
ST	60mm, straight, taper point