


X-Twist Fixation System

Screw-In Suture Anchors

The X-Twist™ is a knotless, open architecture, and fast inserting suture anchor system with multiple sliding suture and suture tape configurations.

The X-Twist system is designed for procedural versatility: single or double row, in medial and/or lateral row applications, and in both knotless and fully knotted constructs, to reproducibly and confidently repair soft tissue to bone.



X-Twist Knotless Screw-In Suture Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material
XT24700	4.75	23	X-Twist PEEK 4.75mm, w/1, #2 suture (blu/wht)	PEEK
XT24701	4.75	23	X-Twist PEEK 4.75mm, w/1, 1.6mm suture tape (blu/wht)	PEEK
XT24702	4.75	23	X-Twist PEEK 4.75mm, w/2, 1.6mm suture tapes (blu, blu/wht)	PEEK
XT24703	4.75	23	X-Twist PEEK 4.75mm, w/2, 1.6mm suture tapes (blk, blk/wht)	PEEK
XT25500	5.5	23	X-Twist PEEK 5.5mm, w/1, #2 suture (blu/wht)	PEEK
XT25501	5.5	23	X-Twist PEEK 5.5mm, w/1, 1.6mm suture tape (blu/wht)	PEEK
XT25502	5.5	23	X-Twist PEEK 5.5mm, w/2, 1.6mm suture tapes (blu, blu/wht)	PEEK
XT25503	5.5	23	X-Twist PEEK 5.5mm, w/2, 1.6mm suture tapes (blk, blk/wht)	PEEK
XT26200	6.25	23	X-Twist PEEK 6.25mm, w/1, #2 suture (blu/wht)	PEEK

X-Twist Instrumentation

Part #	Description	Sterile/Non-sterile	Single-use/Reusable
XT04701	X-Twist 4.75-5.5mm awl	Non-sterile	Reusable
XT04700	X-Twist 4.75mm tap	Non-sterile	Reusable
XT05500	X-Twist 5.5mm tap	Non-sterile	Reusable
XT00001	X-Twist sterilization tray	Non-sterile	Reusable

Suture Tape

Part #	Size (mm)	Description	QTY/ Box
20032S	1.6	Parcus Braid® suture tape (wht/blu)	1
20034S	1.6	Parcus Braid® suture tape (wht/blk)	1
20174S	1.6	Parcus Braid® infinity loop (wht/blu) w/ straight needle	1
20042	2	Parcus Braid® suture tape (wht/blu)	6
20042S	2	Parcus Braid® suture tape (wht/blu)	1
20238	2	Anika Braid suture tape (wht/blu)	12
20238S	2	Anika Braid suture tape (wht/blu)	1
20239	2	Anika Braid suture tape (blu)	12
20239S	2	Anika Braid suture tape (blu)	1
20240	2	Anika Braid suture tape (blk/blu)	12
20240S	2	Anika Braid suture tape (blk/blu)	1

X-Twist mean pullout strength is comparable to the Arthrex SwiveLock C.*

*Data on file.

