

Knotless PEEK CF

Push-in Suture Anchors

Knotless PEEK CF Push-In Suture Anchors are recommended for use in both large and small-joint repairs.

Features & Benefits

Made of carbon fiber-reinforced PEEK-OPTIMA

- Non-absorbable, radiolucent, and MR safe
- Modulus of elasticity closely matching cortical bone

Knotless

- Provides a step-saving alternative to conventional "knotted" suture anchors
- Eliminates "knot stacks" associated with soft tissue irritation



Knotless PEEK CF Push-In Suture Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10312	2.8	10.2	28 Knotless, push-in, suture anchor, w/suture passer	PEEK CF	1
10994	3.5	10	35 Knotless, push-in, suture anchor, w/suture passer	PEEK CF	1
10313	3.5	10	35 Knotless, push-in, suture anchor (implant only, driver not included), w/suture passer	PEEK CF	1
10314	4.5	16.3	45 Knotless, push-in, suture anchor (implant only, driver not included), w/suture passer	PEEK CF	1
11175	4.5	16.3	45 Knotless, push-in, suture anchor, w/suture passer	PEEK CF	1
10315	5.5	18.1	55 Knotless, push-in, suture anchor (implant only, driver not included), w/suture passer	PEEK CF	1
11176	5.5	18.1	55 Knotless, push-in, suture anchor, w/suture passer	PEEK CF	1

Knotless PEEK CF Push-In Suture Anchor, Hip Length

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10312H	2.8	10.2	28 Knotless, push-in, suture anchor, w/hip length driver, w/suture passer	PEEK CF	1

Instrumentation for Knotless PEEK CF Push-In Suture Anchors

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10344	-	Shoulder	35 Knotless suture anchor driver ¹	Sterile	Single-use
10483	-	Shoulder	35 Knotless suture anchor driver ¹	Unsterile	Reusable
10345	-	Shoulder	45 or 55 knotless suture anchor driver ²	Sterile	Single-use
10487	-	Shoulder	45 or 55 knotless suture anchor driver ²	Unsterile	Reusable
10365	3.2	Shoulder	Awl	Unsterile	Reusable
10333U	2.8	Small Joint	Drill bit, w/positive stop	Unsterile	Reusable
10874U	-	Small Joint	28 Knotless drill bit, w/positive stop	Unsterile	Reusable
10538U	2.8	Shoulder	Drill bit, w/positive stop	Unsterile	Reusable
10615U	-	Shoulder	28 Knotless drill bit, w/positive stop	Unsterile	Reusable
10482U	3	Small Joint	Drill bit, w/positive stop	Unsterile	Reusable
10537U	3	Shoulder	Drill bit, w/positive stop	Unsterile	Reusable
10505U	3.2	Small Joint	Drill bit, w/positive stop	Unsterile	Reusable
10540U	3.2	Shoulder	Drill bit, w/positive stop	Unsterile	Reusable
10823	-	Small Joint	8-point slotted drill guide, w/beaver tail handle	Unsterile	Reusable
10924	5.2	Shoulder	8-point slotted knotless drill guide, w/beaver tail handle ³	Unsterile	Reusable
10446	-	Shoulder	8-point drill guide, w/trocar tip obturator	Unsterile	Reusable
10330	-	Shoulder	V-Mouth drill guide w/trocar tip obturator	Unsterile	Reusable
10406	3.4	Shoulder	Blunt tip obturator ⁴	Unsterile	Reusable
10771	3.3	Shoulder	Conical tip obturator, w/1.8mm cannulation ⁴	Unsterile	Reusable
10405	3.5	Shoulder	Trocar tip obturator ⁴	Unsterile	Reusable
10397	2.8	Shoulder	Switching stick	Unsterile	Reusable

Parcus Medical, LLC.

6423 Parkland Dr, Sarasota, FL
 1-941-755-7965 · www.parcusmedical.com
ParcusCS@anika.com

Parcus Medical has joined Anika

